

TLi - WRINKLE REPAIR

This Topical Light Infusion[®] (TLi) Wrinkle Repair program is 21 minutes in duration. Adding **Lightwave Red Light Therapy*** will increase the treatment time. The goal of this program is to drive in powerful ingredients to stimulate a biological change within the skin. Our wrinkle repair serum and cream use a multitude of ingredients including cutting edge peptides, to boost collagen and elastin formation. This potent two step formulation process is designed to reduce flaccidity, lifts sagging tissues, improve skin hydration and restores skin tone and structure resulting in an overall noticeable skin improvement. This program does not include any light therapy benefits.

Note: To activate the Topical Light Infusion[®] setting on the system, you must first purchase Topical Light Infusion[®] topicals. Only Topical Light Infusion[®] topicals can be used with this system. Additionally, take all necessary before pictures of target area prior to starting session number one. If the photos are taken at the conclusion of the treatment, the skin can temporarily appear flushed and red in color.

1. Prepare the Skin: It is important to properly prep the skin for Topical Light Infusion[®]. This is done by thoroughly cleansing the skin with the abi Enzyme Wash. This will ensure all excess oil, make-up, dirt, and debris has been removed. For significantly aged and/or sun damaged skin, additional exfoliation is recommended prior to starting the Topical Light Infusion[®] treatment. Microdermabrasion or abi's various scrubs and chemical peels are all highly effective in properly preparing the skin.

**Lightwave LED Light therapy option: Including an LED light therapy session with your Topical Light Infusion[®] treatment will provide the best possible results. This can be done by completing a standard LIGHTWAVE treatment before completing the Topical Light Infusion[®] program. If you are limited on time, you can pre-treat the skin by turning on the "LW" button at the bottom of your screen prior to starting the TLi program. This will add a RED ONLY LED light session to the beginning of your treatment and increase your treatment time.*

2. Topical Application: Step 1 of 2, remove cap from Step 1 Wrinkle Repair serum syringe and disperse a small amount onto the tips of your index and fore finger. Apply serum on the skin starting at the forehead and working your way down the face to ensure an even amount is applied. Reapply serum to fingertips as needed. You have been provided with enough serum to cover the face and neck area. Allow serum to penetrate the skin for 30 – 60 seconds before moving on to step 2.

Next, Step 2 of 2, remove cap from Step 2 Wrinkle Repair

cream syringe and disperse a small amount onto the tips of your index and for finger. Apply cream on the skin starting at the forehead and working your way down the face to ensure an even amount is applied.

Precautions: Once both topicals have been applied, shield the patient's eyes by fitting the protective goggles in place. The use of goggles is required to minimize brightness and to avoid any incidental eye exposure.

3. Treatment Instructions: Place the LED panel directly over the face and neck area. To maximize the results of the treatment, make sure all six panels are as close to the target area as comfortably possible. Activate the TLi treatment. The LED dome panel will remain lit for 21 minutes. The light is pulsing even though it appears to be constant. Therefore, individuals with epilepsy should be monitored closely.

4. Post Care: Finally, at the end of the Topical Light Infusion treatment, remove the LED panel. Next, apply a dime size amount of the abi Hydrating Recovery Gel to help nourish the skin and dispose of free radicals. Allow the product to penetrate. Then apply a thin layer of ABI Moisturizing Antioxidant to help protect the skin's natural moisture factor. Finally, apply sunscreen to protect the skin from any future UV damage. It is also important to always have the patient increase their water intake by at least 8oz following a treatment.

NUMBER OF SUGGESTED TREATMENTS: An initial eight treatment series is recommended for this protocol over a 4-5 week period. The client should receive Topical Light Infusion[®] treatments two times a week with at least 36 hours between treatment sessions. A follow-up treatment session should occur on week 8-10 which will allow the body's new healthy skin cells to reach the skin's surface and give an accurate reflection of the treatment outcome. Maintenance treatments can prolong the youthful appearance once every 6-8 weeks or as needed. The Combination of LIGHTWAVE™ and Topical Light Infusion[®] are highly recommended for the best possible outcome.